



Lempert Cutter - Angled Down/11cm/9mm Cutting Edge



One edge is finely serrated.

This product is discontinued and has limited availability.

General Info

Product Number	16141-11
Category	Cutters > Medical
EAN	4251909505677



Find 16141-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	9 mm		
Length	11 cm		
Jointaction	Single		
Tipshape	Angled Down		