



Lempert Cutter - Angled Up/11cm/9mm Cutting Edge



One edge is finely serrated.

This product is discontinued and has limited availability.

General Info

Product Number	16142-11
Category	Cutters > Medical
EAN	4251909505684



Find 16142-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cuttingedge	9 mm		
Length	11 cm		
Jointaction	Single		
Tipshape	Angled up		