



Spring Scissors Style Cutter - Straight/Blunt/13cm/12mm Cutting Edge



Cuts fine bone and cartilage.

General Info

Product Number	16144-13
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505707



Find 16144-13 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	12 mm		
Length	13 cm		
Serrated	No		
Tipdiameter	0.4 mm		
Tips	Blunt/Blunt		
Tipshape	Straight		