



Littauer Cutter - Straight/15cm/19mm Cutting Edge



General Info

Product Number	16152-15
Category	Cutters > Medical
EAN	4251909505738



Find 16152-15 on
finescience.com

Technical Info

Shipping Weight	115 g	Net Weight	102 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	19 mm		
Length	15 cm		
Jointaction	Single		
Tipshape	Straight		