



Friedmann-Pearson Rongeur - Curved/0.5mm Cup Diameter/14cm



Exclusive
Line



MOST
POPULAR

Experience the unparalleled precision of our flagship Rongeur from the exclusive FST line. Crafted to perfection, this instrument features a specialized 0.5mm cup tailored for the demanding field of neurological research. Engineered to meet the exacting standards of brain and spinal cord procedures, its razor-sharp cup guarantees effortless cutting while ensuring captured tissue remains securely within the tool. Perfect your procedures with confidence, supported by innovation that sets new benchmarks in surgical excellence. Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.

General Info

Product Number	16221-14
Category	Exclusive Line Rongeurs, Punches & Sampling > Rongeurs Most Popular > Rongeurs, Punches & Sampling
EAN	4251909505752



Find 16221-14 on
finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	53 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	0.5 mm		
Length	14 cm		
Jointaction	Single		
Tipshape	Curved		