




Barraquer Mini Retractor - Wire/2cm/1cm Maximum Spread



General Info

Product Number	17000-01
Category	Retractors Most Popular > Retractors
EAN	4251909505769



Find 17000-01 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	2 cm		
Maximumspread	1 cm		
Teeth	Wire, Blunt		
Teethdepth	4 mm		
Teethlength	3 mm		
Teethwidth	3 mm		