



Goldstein Mini Retractor - 2x2 Blunt/2cm/1cm Maximum Spread



General Info

|                |               |
|----------------|---------------|
| Product Number | 17005-03      |
| Category       | Retractors    |
| EAN            | 4251909505851 |



Find 17005-03 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |     |
|------------------------|---------------------------------|-----------------|-----|
| Shipping Weight        | 20 g                            | Net Weight      | 2 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA  |
| Minimum Order Quantity | 1                               |                 |     |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |     |
| Length                 | 2 cm                            |                 |     |
| Maximumspread          | 1 cm                            |                 |     |
| Teeth                  | 2 x 2 Blunt                     |                 |     |