



Spring Loaded Retractor - Blades/6cm/3cm Maximum Spread



General Info

Product Number	17006-05
Category	Retractors
EAN	4251909505875



Find 17006-05 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	6 cm		
Maximumspread	3 cm		
Teeth	Atraumatic Blades		
Teethdepth	13 mm		
Teethlength	15 mm		
Teethwidth	1 mm		