



Alm Retractor - 4x4 Blunt/7cm/4cm Maximum Spread/16mm Teeth Length



General Info

Product Number	17009-07
Category	Retractors
EAN	4251909505905



Find 17009-07 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	7 cm		
Maximumspread	4 cm		
Shape	Slightly curved down		
Teeth	4 x 4 Blunt		
Teethdepth	5 mm		
Teethlength	3 mm		