

F · S · T

Liebermann Retractor - Blades/8cm/4.5cm Maximum Spread



General Info

Product Number	17009-08
Category	Retractors
EAN	4251909505912



Find 17009-08 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	8 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8 cm		
Feature	Blades		
Maximumspread	4.5 cm		
Shape	Slightly curved down		
Teeth	None		
Teethdepth	4.5 mm		
Teethlength	5.5 mm		