



Heiss Retractor - 4x4 Blunt/10cm/2.3cm Maximum Spread



General Info

Product Number	17011-10
Category	Retractors
EAN	4251909505936



Find 17011-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Maximumspread	2.3 cm		
Teeth	4 x 4 Blunt		
Teethdepth	4 mm		
Teethlength	3 mm		
Teethwidth	3 mm		