



Weitlaner Retractor - 3x4 Blunt/13cm/6.5cm Maximum Spread



General Info

Product Number	17012-13
Category	Retractors
EAN	4251909505950



Find 17012-13 on
finescience.com

Technical Info

Shipping Weight	85 g	Net Weight	81 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Maximumspread	4 cm		
Shape	Straight		
Teeth	3 x 4 Blunt		
Teethlength	16 mm		
Teethwidth	19 mm		