



Wullstein Retractor - 3x3 Blunt/10.5cm/4.5cm Maximum Spread



General Info

Product Number	17018-11
Category	Retractors
EAN	4251909505974



Find 17018-11 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	49 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Maximumspread	3 cm		
Shape	Angled down		
Teeth	3 x 3 Blunt		
Teethlength	16 mm		
Teethwidth	2 mm		