



Fomon Retractor - 2 Pronged Ball/17cm



2mm ball diameter at each tip.

General Info

Product Number	17020-17
Category	Retractors
EAN	4251909505981



Find 17020-17 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	17 cm		
Teeth	2 with 2 mm Ball each tip		
Teethdepth	6 mm		
Teethwidth	13 mm		