



Guthrie Retractor - 2 Pronged Blunt/13cm



General Info

Product Number	17021-13
Category	Retractors
EAN	4251909505998



Find 17021-13 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Teeth	2 Blunt		
Teethdepth	3 mm		
Teethwidth	2.5 mm		