



Blair Retractor - 4 Pronged Blunt/13cm



General Info

Product Number	17022-13
Category	Retractors
EAN	4251909506001



Find 17022-13 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Teeth	4 Blunt		
Teethdepth	3.5 mm		
Teethwidth	6 mm		