



Moria Barraquer Retractor - 3.3 cm/1.5 cm



General Info

Product Number	17300-03
Category	Retractors
EAN	4251909512415



Find 17300-03 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	3.3 cm		
Maximumspread	1.5 cm		
Teeth	Wire		
Teethdepth	5 mm		
Teethlength	10 mm		
Teethwidth	5 mm		