




Moria Barraquer Retractor - 3.5 cm/2 cm



General Info

Product Number	17300-04
Category	Retractors
EAN	4251909512422



Find 17300-04 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	3.5 cm		
Maximumspread	2 cm		
Tapewidth	5 mm		
Teeth	Wire		
Teethdepth	7 mm		
Teethlength	13 mm		