



Small Vessel Cauterizer Kit



This small wireless cauterizer keeps bleeding to a minimum, providing a surgical field clear of clamps and hemostats. Temperature at the tips is approximately 1000C. The unit should be activated only momentarily so the delicate tips to not burn out. Replacement cautery tips are available. Remove batteries before autoclaving. Kit includes: 1 Cautery Handle 1 Straight Fine Tip 1 Angled Fine Tip 2 AA Batteries 1 Storage Case

General Info

Product Number	18000-00
Category	Lab Accessories > Cauterizer
EAN	4251909506032



Find 18000-00 on finescience.com

Technical Info

Shipping Weight	265 g	Net Weight	258 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Handlelength	13.5 cm		
Handleweight	100 g		