



Cannulation System



This micro cannula is designed for quick and easy insertion into arteries and veins. It can also be used as an indwelling catheter during long-term experiments. There is a contoured tip for easy insertion, a shoulder for tie-in to prevent accidental removal, and a female luer connection for linkage to infusion lines. With care it can be reused several times. Detailed instructions are included.

General Info

Product Number	18000-10
Category	Lab Accessories > Cannulas
EAN	4251909506087



Find 18000-10 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) trocar & teflon		
Length	8 cm		
Gauge	27		
Innerdiameter	0.16 mm		
Outerdiameter	0.33 mm		