



Circular Saw for Micro Drill

For cutting bone and other hard substances. Attaches to mandrel (FST No. 18000-18).



General Info

Product Number	18000-20
Category	Lab Accessories > Micro Drill Accessories
EAN	4251909506117



Find 18000-20 on finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Diameter	16 mm		
Shaftdiameter	2.3 mm		
Thickness	0.1 mm		