



Skin Marker

This surgical Skin Marker has a versatile tip for making fine (0.02mm) lines, or a broad (2mm) line when held at an angle. The ink is violet, non toxic, non tattooing and surgical prep resistant. Markings can be removed from the skin with isopropyl (rubbing) alcohol. A six-inch ruler calibrated in both centimeters and inches is included.



General Info

Product Number	18000-30
Category	Lab Accessories > Marking & Identification
EAN	1088776100500



Find 18000-30 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Various		
Color	violet		