

Used for precise measurement during surgery.



Allgemein

Produktnummer	18000-35
Kategorie	Lab Accessories > Measurement
EAN	4251909506124



Finde 18000-35 auf
finescience.com

Technische Informationen

Versandgewicht	50 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8.5 cm		
Increments	0.5 mm		
Maximummeasurement	20 mm		
Tipdiameter	0.5 mm		