



Castroviejo Surgical Caliper

Used for precise measurement during surgery.



General Info

Product Number	18000-35
Category	Lab Accessories > Measurement
EAN	4251909506124



Find 18000-35 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8.5 cm		
Increments	0.5 mm		
Maximummeasurement	20 mm		
Tipdiameter	0.5 mm		