



Trephines for Micro Drill - 2.7mm Tip Diameter/Stainless Steel



These trephines with a side opening are ideal for clearing bone and other materials.

General Info

Product Number	18004-27
Category	Lab Accessories > Micro Drill Accessories
EAN	4251909506186



Find 18004-27 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	44 mm		
Shaftdiameter	2.3 mm		
Tipdiameter	2.7 mm		