



This high temperature, heat cautery is a great alternative to disposable cautery units. Temperature at the tips is approximately 1000° C.

- Kit includes:
- 1 Cautery Handle
 - 2 Standard Tips (18010-01)
 - 2 Fine Tips (18010-02)
 - 2 AA Batteries
 - 1 Storage Case

General Info

Product Number	18010-00
Category	Lab Accessories > Cauterizer Most Popular > Lab Accessories
EAN	4251909506209



Find 18010-00 on
finescience.com

Technical Info

Shipping Weight	260 g	Net Weight	258 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Various		
Autoclavesafe	No		
Handlelength	16 cm		
Handleweight	60 g		