



Replacement Tips for High Temp. Cautery Unit, Fine tips



General Info

Product Number	18010-02
Category	Lab Accessories > Cauterizer
EAN	4251909506223



Find 18010-02 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	15 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Sterile	Yes		
Tipdiameter	0.25 mm		
Tipshape	Straight / Fine		