



Replacement Tips for Low Temp. Cautery Unit - Elongated Straight / Fine



Allgemein

Produktnummer	18019-02
Kategorie	Lab Accessories > Cauterizer
EAN	4251909510428



Finde 18019-02 auf finescience.com

Technische Informationen

Versandgewicht	172 g	Nettogewicht	170 g
Paketgröße	10 Pc.	Maßeinheit	PA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Sterile	Yes		
Tipdiameter	0.25 mm		
Tipshape	Elongated Straight / Fine		