



Suture Thread - Braided Silk/Size 1



General Info

Product Number	18020-01
Category	Needles, Threads & Holders > Thread Most Popular > Needles, Threats & Holders
EAN	4251909506247



Find 18020-01 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Braided Silk		
Autoclavesafe	Yes		
Size	1		
Spoollength	23 m		
Sterile	No		
Suturediameter	0.45 mm		
Suturematerial	Braided Silk		