




Suture Thread - Braided Silk/Size 3/0



Allgemein

Produktnummer	18020-30
Kategorie	Needles, Threads & Holders > Thread Most Popular > Needles, Threats & Holders
EAN	4251909506285



Technische Informationen

Versandgewicht	30 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Braided Silk		
Size	3/0		
Spoollength	23 m		
Sterile	No		
Suturediameter	0.22 mm		
Suturematerial	Braided Silk		