




Suture Thread - Braided Silk/Size 7/0



Allgemein

Produktnummer	18020-70
Kategorie	Needles, Threads & Holders > Thread Most Popular > Needles, Threats & Holders
EAN	4251909506322



Technische Informationen

Versandgewicht	20 g	Nettogewicht	16 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Braided Silk		
Autoclavesafe	Yes		
Size	7/0		
Spoollength	91 m		
Sterile	No		
Suturediameter	0.05 mm		
Suturematerial	Braided Silk		