



Tuebingen Suture Tying Forceps - 10cm



A pin inside the shanks limits the applied pressure.

General Info

Product Number	18025-10
Category	Forceps & Tweezers > Smooth
EAN	4251909506339



Find 18025-10 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.5 x 0.5 mm		
Tips	Smooth with platform		
Tipshape	Straight		
Tyingplatform	5 mm		