



Suture Nylon Monofilament Suture - Nylon monofilament/Size 3/0



Nylon is known for its uniform diameter flexibility and outstanding consistent knot security.

General Info

Product Number	18030-30
Category	Needles, Threads & Holders > Thread
EAN	4251909506353



Find 18030-30 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Nylon Monofilament		
Size	3/0		
Spoollength	91 m		
Sterile	No		
Suturediameter	0.2 mm		
Suturematerial	Nylon monofilament		