



Suture Nylon Monofilament Suture - Nylon monofilament/Size 4/0



Nylon is known for its uniform diameter flexibility and outstanding consistent knot security.

General Info

Product Number	18030-40
Category	Needles, Threads & Holders > Thread
EAN	4251909506360



Find 18030-40 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Nylon Monofilament		
Autoclavesafe	Yes		
Size	4/0		
Spoollength	91 m		
Sterile	No		
Suturediameter	0.15 mm		
Suturematerial	Nylon monofilament		