



Sample Corer - 4G



This specialized instrument has been developed for exact quantitative analyses and comes complete with plunger. The precise construction ensures reproducible results and reliable measurement accuracy in laboratory analytical applications.

General Info

Product Number	18035-05
Category	Rongeurs, Punches & Sampling > Sampling
EAN	4251909506445



Find 18035-05 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	33 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8.5 cm		
Gauge	4		
Innerdiameter	5 mm		
Outerdiameter	6 mm		