



Sample Corer - 17G



This specialized instrument has been developed for exact quantitative analyses and comes complete with plunger. The precise construction ensures reproducible results and reliable measurement accuracy in laboratory analytical applications.

General Info

Product Number	18035-80
Category	Rongeurs, Punches & Sampling > Sampling
EAN	4251909506476



Find 18035-80 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8.5 cm		
Gauge	17		
Innerdiameter	0.8 mm		
Outerdiameter	1.45 mm		