



For removal of frozen brain microregions.

### General Info

Product Number	18036-19
Category	Rongeurs, Punches & Sampling > Punches
EAN	4251909506483



Find 18036-19 on  
[finescience.com](https://www.finescience.com)

### Technical Info

Shipping Weight	1 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	4.5 cm		
Gauge	19		
Innerdiameter	0.69 mm		
Outerdiameter	1 mm		