



NIH Style Neuro Punch - 24G



For removal of frozen brain microregions.

General Info

Product Number	18036-24
Category	Rongeurs, Punches & Sampling > Punches
EAN	4251909506506



Find 18036-24 on
finescience.com

Technical Info

Shipping Weight	1 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	4.2 cm		
Gauge	24		
Innerdiameter	0.31 mm		
Outerdiameter	0.61 mm		