



Provides good access to deep seated vessels.

General Info

Product Number	18038-45
Category	Clamps & Clips > Clamps
EAN	4251909506513



Find 18038-45 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	8 mm		
Clamppress	Light		
Length	52.5 mm		
Grampress	80 g		
Jawdims	8.6 x 1.3 mm		
Serrated	Yes		
Shape	Straight		