




Double Micro Clamp - Straight/11mm

For vascular approximation at anastomoses. These clamps are mounted on a 30 mm sliding bar and can easily be moved by applying gentle force. The distance between the jaws of the clamps is adjustable from 3 to 16 mm.



General Info

Product Number	18040-22
Category	Clamps & Clips > Clamps
EAN	4251909506568



Find 18040-22 on finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5.5 mm		
Clamppress	Variable		
Length	11 mm		
Grampress	5 to 15 g		
Jawdims	5.5 x 1.5 mm		
Serrated	No		
Shape	Straight		