



Dieffenbach Clamp - Straight/50mm



MOST
POPULAR



SERRATED

General Info

Product Number	18050-50
Category	Clamps & Clips > Serrefines Most Popular > Clamps & Clips
EAN	4251909506605



Find 18050-50 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15 mm		
Clamppress	Strong		
Length	50 mm		
Jawdims	15 x 2 mm		
Serrated	Yes		
Shape	Straight		