



Dieffenbach Clamp - Curved/32mm



NEW ITEM



SERRATED

General Info

Product Number	18051-32
Category	Clamps & Clips > Serrefines
EAN	4251909512590



Find 18051-32 on
finescience.com

Technical Info

Shipping Weight	3 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Clampinglength	10 mm		
Clamppress	Medium		
Length	32 mm		
Jawdims	10 x 2 mm		
Serrated	Yes		
Shape	Curved		