



Dieffenbach Clamp - Curved/50mm



Allgemein

Produktnummer	18051-50
Kategorie	Clamps & Clips > Serrefines Most Popular > Clamps & Clips
EAN	4251909506636



Finde 18051-50 auf
finescience.com

Technische Informationen

Versandgewicht	5 g	Nettogewicht	2 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15 mm		
Clamppress	Strong		
Length	50 mm		
Jawdims	15 x 2 mm		
Serrated	Yes		
Shape	Curved		