



V-Hook - 11cm/0.4mm Tip Diameter



Aids in feeding fine suture around blood vessels, works very similar to ligation aids.

General Info

Product Number	18052-12
Category	Needles, Threads & Holders
EAN	4251909506674



Find 18052-12 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Holediameter	1 mm		
Tipdimensions	0.4 x 0.2 mm		