



Smallest bulldog type hemostatic clamps with a 3mm spring width.

General Info

Product Number	18053-28
Category	Clamps & Clips > Serrefines
EAN	4251909506681



Find 18053-28 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Clamping length	6 mm		
Clamp press	Medium		
Length	28 mm		
Jaw dims	6 x 2 mm		
Serrated	Yes		
Shape	Curved		