



Mini Serrefine - Straight/28mm



Smallest bulldog type hemostatic clamps with a 3mm spring width.

General Info

Product Number	18054-28
Category	Clamps & Clips > Serrefines
EAN	4251909506698



Find 18054-28 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	6 mm		
Clamppress	Medium		
Length	28 mm		
Jawdims	6 x 2 mm		
Serrated	Yes		
Shape	Straight		