




Micro-serrefine with delicate, atraumatic serrations. Used to obtain a clear surgical field for extended arterial cross-clamping. Caution: These clips should never be compressed with fingers. Use only the appropriate FST applicator 18056-14 or 18057-14.

General Info

Product Number	18055-03
Category	Clamps & Clips > Serrefines Most Popular > Clamps & Clips
EAN	4251909506728



Find 18055-03 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	6 mm		
Clamppress	Medium		
Length	15 mm		
Grampress	100 g		
Jawdims	6 x 1 mm		
Serrated	Yes		
Shape	Straight		