




Micro-serrefine with delicate, atraumatic serrations. Used to obtain a clear surgical field for extended arterial cross-clamping. Caution: These clips should never be compressed with fingers. Use only the appropriate FST applicator 18056-14 or 18057-14.

General Info

Product Number	18055-04
Category	Clamps & Clips > Serrefines Most Popular > Clamps & Clips
EAN	4251909506735



Find 18055-04 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	4 mm		
Clamppress	Strong		
Length	13 mm		
Grampress	125 g		
Jawdims	4 x 0.75 mm		
Serrated	Yes		
Shape	Straight		