



Biemer Clip - 6 mm/15 mm



Allgemein

Produktnummer	18058-15
Kategorie	Clamps & Clips > Clamps MedTech Line
EAN	4251909511562



Finde 18058-15 auf
finescience.com

Technische Informationen

Versandgewicht	5 g	Nettogewicht	1 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	6 mm		
Clamppress	Light		
Length	15 mm		
Grampress	0.29 - 0.39 N / 30 - 40 g		
Serrated	Yes		
Shape	Straight		