



Allgemein

Produktnummer	18060-18
Kategorie	Lab Accessories > Feeding Needles Most Popular > Lab Accessories
EAN	4251909511807



Finde 18060-18 auf
finescience.com

Technische Informationen

Versandgewicht	3 g	Nettogewicht	2 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	38 mm		
Gauge	18		
Shape	Straight		
Sterile	No		
Style	Reusable		
Tipdiameter	1.7 mm		
Tipshape	Barrel		
Xsizeweightmc	20-30 g		
Xsizeweightrt	100-200 g		
Xsizeweightgp	250-300 g		